

THIS FILE IS MADE AVAILABLE THROUGH THE DECLASSIFICATION EFFORTS AND RESEARCH OF:

# THE BLACK VAULT

THE BLACK VAULT IS THE LARGEST ONLINE FREEDOM OF INFORMATION ACT / GOVERNMENT RECORD CLEARING HOUSE IN THE WORLD. THE RESEARCH EFFORTS HERE ARE RESPONSIBLE FOR THE DECLASSIFICATION OF THOUSANDS OF DOCUMENTS THROUGHOUT THE U.S. GOVERNMENT, AND ALL CAN BE DOWNLOADED BY VISITING:

[HTTP://WWW.BLACKVAULT.COM](http://www.blackvault.com)

YOU ARE ENCOURAGED TO FORWARD THIS DOCUMENT TO YOUR FRIENDS, BUT PLEASE KEEP THIS IDENTIFYING IMAGE AT THE TOP OF THE .PDF SO OTHERS CAN DOWNLOAD MORE!



REPLY TO  
ATTENTION OF

DEPARTMENT OF THE ARMY  
U.S. ARMY CHEMICAL AND BIOLOGICAL DEFENSE COMMAND  
5232 FLEMING ROAD  
ABERDEEN PROVING GROUND, MARYLAND 21010-5423

August 25, 1998

Freedom of Information and Privacy Act Office

Mr. John Greenewald, Jr.

Dear Mr. Greenewald:

I am writing in response to your June 4, 1998, Freedom of Information Act (FOIA) request.

Enclosed you will find a copy of "The Determination of Trace Quantities of BZ (3-Quinuclidinyl Benzilate) In Air and Water".

We have incurred the following reimbursable fees:

Professional Hours (1/2 hr at \$25.00 per hour) (searching)	\$12.50
--	---------

Office Copy Reproduction (56 pages at 15 cents per page)	\$ 8.40
---	---------

**TOTAL FEES DUE: \$20.90**

Please forward your check to me in the amount of \$20.90, payable to Treasurer of the United States, at the following address:

Mrs. Cheryl S. Fields  
US Army Chemical & Biological Defense Command  
Attention: AMSCB-CIH, Building E5232  
Aberdeen Proving Ground, Maryland 21010-5423

As a category iii requester (other), you are entitled to 100 free pages. However, you have already received more than 100 pages of documents free; therefore, fees are payable upon receipt of this letter. The enclosed document is less than 100 pages, but all documents that you receive pertaining to BZ will be compiled and considered as an aggregated request.

I look forward to processing other FOIA requests for you. If you have any questions, please feel free to contact me at (410) 436-1288.

Sincerely,

A handwritten signature in cursive script, reading "Cheryl S. Fields". The signature is written in dark ink and is positioned above the printed name.

Cheryl S. Fields

Freedom of Information and  
Privacy Act Officer

Enclosure

#### Disclaimer

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

#### Disposition

Destroy this report when no longer needed. Do not return it to the originator.


UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

ARCSL-TR-83031

## REPORT DOCUMENTATION PAGE

\*ADB 079073 002\*

1. REPORT NUMBER ARCSL-TR-83031		2. 	
4. TITLE (and Subtitle) The Determination of Trace Quantities of BZ (3-Quinuclidinyl Benzilate) in Air and Water		5. TYPE OF REPORT & PERIOD COVERED Technical Report April 1978 - February 1980	
		6. PERFORMING ORG. REPORT NUMBER	
7. AUTHOR(s) Ellen S. Vigus Albert Deiner		8. CONTRACT OR GRANT NUMBER(s)	
9. PERFORMING ORGANIZATION NAME AND ADDRESS Commander, Chemical Systems Laboratory ATTN: DRDAR-CLN-D Aberdeen Proving Ground, Maryland 21010		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS Project C23	
11. CONTROLLING OFFICE NAME AND ADDRESS Commander, Chemical Systems Laboratory ATTN: DRDAR-CLJ-IR Aberdeen Proving Ground, Maryland 21010		12. REPORT DATE June 1983	
		13. NUMBER OF PAGES 55	
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		15. SECURITY CLASS. (of this report) UNCLASSIFIED	
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE NA	
16. DISTRIBUTION STATEMENT (of this Report) Distribution limited to US Government agencies only because of test and evaluation; June 1983. Other requests for this document must be referred to Commander, Chemical Systems Laboratory, ATTN: DRDAR-CLJ-IR, Aberdeen Proving Ground, Maryland 21010.			
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)			
18. SUPPLEMENTARY NOTES			
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) BZ 3-Quinuclidinyl Benzilate Incapacitating Agents Demilitarization Gas Chromatography Aerosol Particulates Hydrolysis Analysis			
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) Methods were developed for the analysis of BZ (3-quinuclidinyl benzilate) in the emissions and waste materials from demilitarization processes using gas chromatography.  BZ was quantitatively removed from aerosol filters using acidic aqueous solvents. Extraction of BZ from water samples and from the aerosol filter ex- tracts involved the use of chlorinated solvents. Recoveries of BZ at different concentrations and in the presence of possible interferents are given.			

## PREFACE

The work described in this report was authorized under project C23, Disposal of Incapacitating BZ/Agent Munitions and funded by US Army Toxic and Hazardous Materials Agency. The work was started in April 1978 and was completed in February 1980. The experimental data are recorded in Chemical Systems Laboratory notebooks 9787, 9811, and 9926.

The use of trade names in this report does not constitute an official indorsement or approval of the use of such commercial hardware or software. This report may not be cited for purposes of advertisement.

Reproduction of this document in whole or in part is prohibited except with permission of the Commander, Chemical Systems Laboratory, ATTN: DRDAR-CLJ-IR, Aberdeen Proving Ground, Maryland 21010. However, the Defense Technical Information Center is authorized to reproduce the document for United States Government purposes.

This report has not been approved for release to the public.

### Acknowledgement

The authors acknowledge the assistance of Mr. Raymond E. Herd, Research Division, Chemical Systems Laboratory, for the preparation of BZ standard solutions.

Blank

# CONTENTS

	Page
1. INTRODUCTION .....	7
2. MATERIALS AND METHODS .....	9
2.1 Aerosol Studies .....	9
2.1.1 Equipment .....	9
2.1.2 Gas Chromatographic Conditions .....	9
2.1.3 Reagents .....	10
2.1.4 Purification Procedures .....	10
2.2 Water Studies .....	10
2.2.1 Equipment .....	10
2.2.2 Gas Chromatographic Conditions .....	11
2.2.3 Reagents .....	11
2.2.4 Purification Procedures .....	12
3. INVESTIGATIONAL PROCEDURES .....	13
3.1 Aerosol Studies .....	13
3.1.1 Preparation of BZ Standards .....	13
3.1.2 Internal Standard Interferences .....	13
3.1.3 Extraction of BZ from Filters .....	13
3.1.4 Aeration Study .....	14
3.1.5 Concentration Study .....	15
3.2 Water Studies .....	15
3.2.1 Preparation of BZ Standards .....	15
3.2.2 Preparation of Michler's Ketone Standards .....	15
3.2.3 Extraction of BZ from Tap Water .....	15
3.2.4 Extraction of BZ from Tap Water in the Presence of Contaminants .....	17
4. RESULTS .....	17
4.1 Aerosol Studies .....	17
4.2 Water Studies .....	19
5. DISCUSSION OF RESULTS .....	37
5.1 Aerosol Studies .....	37
5.1.1 Statistical Analysis of the Data .....	37
5.1.1.1 Peak Height Measurements .....	37
5.1.1.2 Sample Size .....	40
5.1.1.3 Area Measurement .....	41
5.1.1.4 t-Test and F-Test .....	41
5.2 Water Studies .....	43
6. CONCLUSIONS .....	47
6.1 Employment of NPD for BZ Analysis .....	47
6.2 Employment of FID for BZ Analysis .....	48
LITERATURE CITED .....	49
DISTRIBUTION LIST .....	51



Blank

THE DETERMINATION OF TRACE QUANTITIES OF BZ  
(3-QUINUCLIDINYL BENZILATE) IN AIR AND WATER

1. INTRODUCTION

BZ (3-quinuclidinyl benzilate) was introduced into the United States chemical agent munitions inventory in the early 1960's. In late 1976, a program was initiated by the US Army Toxic and Hazardous Materials Agency for the demilitarization of the remaining stores of bulk agent, the M43 and M44 munitions, contaminated residues, and scrap.

The most promising disposal methods outlined in the Concept Plan<sup>1</sup>, prepared for the demilitarization action, were incineration and chemical neutralization (hydrolysis or transesterification). Since the ultimate disposal method was not specified, a sensitive and specific analytical method was needed that could analyze for BZ in the emissions and waste materials produced by either disposal method.

Based on earlier work on the analysis of BZ, the most promising methods appeared to be gas chromatography (GC) or spectroscopy. GC had been investigated by both Sass<sup>2</sup> and Ellin<sup>3</sup>. Sass demonstrated that GC could be used to separate and quantitate BZ in the presence of known process intermediates. Ellin was able to detect BZ in such complex media as blood and urine by a simple extraction and concentration technique. However, the levels analyzed in both studies were at least 1000 times greater than the levels required by current emission standards. Spectroscopic methods such as Tropaeolin 00<sup>4</sup> and iodine complexation<sup>5</sup> were shown to have sufficient sensitivity, but suffered from interference from the products of hydrolysis (e.g., 3-quinuclidinol).

To determine which of the available methods should be developed further, an evaluation was performed by Battelle Columbus Laboratories (BCL) under contract number DAAK40-73-C-0142<sup>6</sup>. This evaluation concluded that the best choices would be GC or gas chromatography/mass spectrometry (GC/MS). It was found that GC, equipped with a flame ionization, nitrogen/phosphorus, or electron capture detector had sufficient sensitivity for the analysis of BZ at the emission levels. At the time, it was anticipated that GC could be used for stack, waste stream, and wash-down samples after a simple extraction. It was also felt that a simple cleanup procedure would be required for more complex samples such as hydrolysis brines. It was postulated that GC/MS, being more specific, could be used to detect BZ even in complex brine samples after a simple extraction into an organic solvent.

From the BCL study, GC was chosen over the spectroscopic techniques as the method to be developed for most of the routine analyses needed in the BZ demilitarization program. The nitrogen-phosphorus detector (NPD) offered the most selectivity of the detectors evaluated by BCL for routine analyses. Detector selectivity was desirable in that it might eliminate the need for any complex cleanup methods. The GC/MS was to be used only for very complex samples that could not be analyzed by simpler GC detectors and as a referee method for more specific analysis of BZ in any questionable samples.

As part of the BCL evaluation, the question arose as to whether BZ would be lost from aerosol filters during sampling. In one test, less than 10 percent of 1  $\mu$ g of a BZ spike applied to a filter was recovered after aerating with 500 liters of air. Losses of this magnitude would have seriously hampered the development of analytical methods for BZ particulates.

Initiated shortly after contract DAAK40-73-C-0142 was completed, a study at Chemical Systems Laboratory (CSL) first looked at aeration losses from filters. Based on the BCL methods evaluated, GC with an NPD was chosen as the analytical method for this work. The BZ levels and aeration times used in the CSL work<sup>7</sup> were based on the then existing chronic exposure monitoring requirements, 4 mg-min/ $m^3$ . Assuming an 8-hour workday, the permissible chronic exposure concentration would be  $8.3 \times 10^{-3}$  mg/ $m^3$ . If a real-time monitor (RTM) were assumed to operate at 10 l/min for 10 minutes, the amount of BZ on the filter at the maximum permissible concentration would be 0.8  $\mu$ g.

No firm requirements were available at the time for a personnel monitor. It was assumed that a level of one-tenth the requirement for a RTM would approximate the level for this purpose. The levels used in the filter study were set at the maximum levels for the RTM and the personnel monitor. The aeration volumes were set at the assumed total volumes for the RTM and the personnel monitor.

In 1978, BCL was awarded a second contract (DAAK11-78-C-0096) for conducting the laboratory phase studies required for the design of the BZ munitions demilitarization facility. Part of this work involved the development and verification of analytical methods to support BZ agent incineration and neutralization tests. BCL was also tasked to prepare analytical methods manuals for eventual plant operation. Early in this work, incineration became the preferred agent destruction process. Analytical tasks were quickly shifted to support an operating plant utilizing incineration equipment. Key analytical methods involved determinations of BZ in furnace effluent gases, as an aerosol particulate in the workplace and in plant washdown/scrubber liquids.

Work was performed at CSL to supplement BCL's analytical methods development. CSL's efforts were directed toward developing an analysis for BZ in water.

Problems with the NPD during the filter study at CSL and later work at BCL during a detector evaluation<sup>8</sup> led to the determination that the NPD was too unstable and fragile for routine use. As a result, GC with a flame ionization detector (FID) was chosen for subsequent methods development. The FID was chosen because of its ruggedness, reliability, ease of operation, and wide linear dynamic range.

The BZ levels used in this part of the study were set at one-quarter the maximum permissible BZ concentration (1 ppb), the maximum concentration (4 ppb), and 10 times the maximum concentration (40 ppb)<sup>7</sup>

in water. At the time, it was believed that the aqueous analysis would be a minor part of the analytical program. As the BCL study progressed, however, it became apparent that the methods for air and water would merge into a singularity, namely, extraction from an aqueous matrix into an organic solvent with final quantification of the BZ by GC or by GC/MS.

The final results of the methods development for air samples using the NPD and for water samples using the FID are documented herein.

## 2. MATERIALS AND METHODS

### 2.1 Aerosol Studies.

#### 2.1.1 Equipment:

- Gas chromatograph - Hewlett-Packard Model 5840A with a Hewlett-Packard Model 18847A/8A nitrogen-phosphorus thermionic flame ionization and a 4 ft x 1/8-inch ID glass column packed with 2 percent OV101 on 100/120 Chromsorb WHP®. (The 2 percent OV101 column was used because the NPD needed a low bleed column.)
- Reflux timer - Maximum setting 120 seconds (Eagle Signal Company, Davenport, Iowa).
- Solenoid - Valve number V52DB2052, 10 watts, 50 psi, 3/16-inch orifice (Skinner Electric Valves Division, New Britain, Connecticut).
- Filters - Gelman Type A/E glass fiber.
- Filter holders - Gelman product no. 1220 or the equivalent.
- Syringes - (10 µl, 25 µl, 50 µl, 100 µl) Hamilton Models 701N, 702LT, 705LT, and 710LT.
- Volumetric glassware to conform to NBS Class A specifications.

#### 2.1.2 Gas Chromatographic Conditions:

- Temperatures - Oven - 230°C  
Detector and injection port - 260°C
- Helium - 30 ± 1 ml/min
- Air - 50 ml/min
- Hydrogen - 3 ml/min

### 2.1.3 Reagents:

- BZ - Obtained from Analytical Branch, Research Division. Analysis by a standard titrimetric method gave a purity of 99+ percent.
- 1,4-dioxane - Distilled-in-glass (Burdick and Jackson Laboratories, Muskegon, Michigan) or freshly purified and distilled ACS reagent grade.
- Methanol - ACS reagent grade.
- Isopropanol - ACS reagent grade.
- Sodium hydroxide - ACS reagent grade.
- Toluene - ACS reagent grade.

All other reagents were ACS or CP in quality.

### 2.1.4 Purification Procedures.

If distilled-in-glass, 1,4-dioxane is not available, ACS reagent grade 1,4-dioxane may be purified and distilled by the following procedures developed during the course of this study; this procedure will produce dioxane that is suitable as a chromatographic solvent.

Three liters of 1,4-dioxane were gently refluxed for 6 to 8 hours with 50g of sodium hydroxide pellets, stoppered, and allowed to cool overnight.

The next morning, the material was placed into a distillation apparatus having an 18-inch insulated Vigreux column. The material was brought to boiling and allowed to reflux for 1 hour. Distillation was then started.

All material boiling below 100°C was discarded (approximately 400 ml). The fraction boiling between 100°C and 102°C was collected (approximately 2 liters). The remaining material in the still was discarded.

The collected fraction was placed into pint size, narrow neck, brown glass bottles which were tightly sealed, and frozen in a refrigerator until use.

## 2.2 Water Studies.

### 2.2.1 Equipment:

- Gas chromatograph - Hewlett-Packard Model 5840A with a Hewlett-Packard Model 18812B FID and a 4 ft x 1/8-inch ID glass column packed with 3 percent OV17 on 80/100 Chromsorb WHP®. (The 3 percent OV17 column gave better chromatographic separation than the 2 percent OV101 column. The higher bleed column was not a problem with the FID).





- Toluene - ACS reagent grade.
- Cyclohexane - ACS reagent grade.
- Sodium sulfate - ACS reagent grade.
- Sodium metabisulfite ( $\text{Na}_2\text{S}_2\text{O}_5$ ) - ACS reagent grade.
- Water (double distilled) - distilled  $\text{H}_2\text{O}$  was redistilled from a well-leached, all-glass distillation apparatus.
- 4,4-bis (dimethylamino) benzophenone [Michler's Ketone] - purified and recrystallized as described in 2.2.4.

All other reagents were ACS or CP in quality.

#### 2.2.4 Purification Procedures.

Chloroform - The ethanol which is used as a preservative in chloroform could destroy the BZ by causing it to transesterify. Acids which form as chloroform ages can also react with the BZ. To remove the ethanol and acids, the chloroform was washed with a sodium carbonate solution, then with double-distilled water as follows:

Approximately 1500 ml of chloroform was washed by vigorously stirring for at least 1 hour with 200 ml of 0.5 M aqueous sodium carbonate. The chloroform was then separated and vigorously stirred three times for 10 minutes each with double-distilled water. The fourth and final double-distilled water wash was performed by shaking for 5 minutes in a 2-liter separatory funnel. After the aqueous and chloroform layers separated, the chloroform was drawn off into a brown-glass bottle.

The Michler's Ketone was dissolved in acetone, treated with decolorizing charcoal, and gently boiled for 5 minutes. The solution was filtered while hot to remove the charcoal. The filtrate was heated to boiling and distilled water was added dropwise to the point just before crystallization occurred. The filtrate was then allowed to cool undisturbed. The resulting crystals were vacuum filtered. The recrystallization from acetone and water was repeated once again. The product was vacuum filtered and dried in a vacuum dessicator set at  $40^\circ\text{C}$  for 48 hours. The resulting product was pure white leaflets with a melting point of  $172^\circ\text{C}$ . Merck Index cites  $172^\circ\text{C}$ .<sup>9</sup>

### 3. INVESTIGATIONAL PROCEDURES

#### 3.1 Aerosol Studies.

##### 3.1.1 Preparation of BZ Standards.

BZ stock solution was prepared by dissolving BZ in toluene. The solution was sealed in serum vials with Teflon® lined septa and stored under refrigeration.

The BZ standards used for the analysis by GC were prepared by dilution of the stock solution with dioxane. Dioxane was chosen because BZ was soluble in it, and the dioxane gave a low background on the NPD. The standards were sealed in serum vials with Teflon® lined septa. The standards were prepared fresh daily for the low concentration, and every other day for the high concentration.

##### 3.1.2 Internal Standard Interferences.

A problem was encountered with the internal standard chosen (1-(N-piperidyl) isopropyl benzilate). An impurity in the concentrated extracts had the same retention time as the internal standard on the column used. This made quantification of the BZ using the internal standard impossible. Time limitations precluded a search for another standard, so efforts were made at eliminating the interference.

The Gelman glass fiber type A/E filters were extracted with dioxane for 8 hours in a Soxhlet extractor. The dioxane was discarded and the filters were extracted for 8 hours with isopropanol. The filters were then dried for 16 hours in a heated vacuum desiccator. This extraction eliminated some of the contamination observed in the chromatogram and greatly improved the baseline stability but did not completely remove the impurity interfering with the internal standard.

Both the distilled-in-glass and freshly-purified and distilled dioxane were concentrated in the same manner as the filter extracts to determine if the impurity was in the dioxane. The results did not conclusively point to the dioxane as the source of contamination.

An effort was made to identify the interfering gas chromatographic peak. A concentrated filter extract sample containing this interference was submitted to BCL for gas chromatographic and mass spectrophotometric analysis under contract DAAK11-78-C-009610. Computer identification of the separated interference pointed to a structure similar to trioctylphosphate. Compounds of this type are widely used as plasticizers in spinning fibers and are common pollutants.

The source of contamination was not determined.

##### 3.1.3 Extraction of BZ from Filters.

The first part of the study determined the efficiency of the procedure used for extracting BZ from glass-fiber filters.



For this study, the filters were impregnated with BZ by administering BZ solutions dropwise with a syringe and air drying for 15 minutes. The filters were transferred to a centrifuge tube and then macerated with about 5 ml of dioxane. The suspension was centrifuged for 3 minutes. The supernatant liquid was decanted into a concentration tube. This extraction, maceration, and centrifugation procedure was repeated twice more. The dioxane extracts were combined and concentrated to about 0.5 ml by heating in a constant-temperature water bath set at 80°C while purging with nitrogen. The nitrogen flow was controlled by a reflux timer connected to a solenoid. This alternately started and stopped the nitrogen at 30-second intervals. The purpose of the interrupted flow was to set up a reflux to help wash down the walls of the concentration tube. When the concentrates reached a level between 0.2 ml and 0.5 ml, they were removed from the bath and diluted to a total volume of exactly 0.5 ml with dioxane. The concentrates were stirred on a Vortex mixer before analyzing by GC.

It was found that the syringe used for injecting samples into the GC had to be thoroughly rinsed with methanol between injections to remove the BZ. It was found that injections had to be made immediately after taking a sample to prevent adsorption of the BZ onto the walls of the syringe barrel. It was also found that BZ could be chromatographed without apparent breakdown at temperatures above its decomposition point if the BZ was injected directly onto the column packing and not onto the glass wool or the dead space above the glass wool.

The interference problem precluded quantification of BZ using an internal standard, so the amount of BZ in the concentrates was determined by direct comparison of concentrate peak heights and areas to standard peak heights and areas.

The extraction efficiency was calculated by dividing the amount of BZ recovered by the amount deposited on the filter.

#### 3.1.4 Aeration Study.

The second part of the study compared the recovery of BZ from aerated filters to the recovery from unaerated filters. Two aeration volumes were studied. The low volume (100 liters) represented sampling at a low rate for a long period such as in the case of a personnel monitor. The high volume (1000 liters) represented sampling at a very high rate such as would be done in stack sampling or area monitoring. The low volume filters were aerated at 4.4 l/min for 25 minutes. The high volume filters were aerated at 13.5 l/min for 75 minutes.

For the aeration study, the filters were spiked with BZ and air dried for 15 minutes. The filters were then aerated by pulling room air through them with a vacuum pump. After aeration, the filters were extracted, and the extracts were concentrated and analyzed in the same manner as in the extraction study.

The recovery from the aerated filters was compared to the recovery for the unaerated filters to determine if any loss of BZ occurred during aeration.

### 3.1.5 Concentration Study.

BZ solutions were concentrated in the same manner as the filter extracts to determine if any BZ was lost during the concentration step. The two levels of BZ studied (0.1 µg and 1.4 µg) were the same as in the extraction and aeration studies.

The BZ solutions were prepared by adding 5 ml of dioxane and a known amount of BZ stock solution to a concentration tube. Additional dioxane was then added to bring the final volume to 12 ml. This was done to approximate the final volume of the filter extracts. The solutions were stirred on a Vortex mixer before they were concentrated. The recoveries were calculated by direct comparison of the concentrate peak heights.

### 3.2 Water Studies.

#### 3.2.1 Preparation of BZ Standards.

- BZ (base) - BZ stock solutions were prepared by dissolving BZ in toluene. The solutions were sealed in serum vials using Teflon® lined septa and stored under refrigeration. The BZ standards used for the analysis by GC were prepared by dilution of the stock solution with methyl ethyl ketone. The standards were dissolved in methyl ethyl ketone because it was a readily available solvent in which both the BZ and the Michler's ketone were soluble and stable. The standards were sealed in serum vials using Teflon® lined septa.

- BZ·HCl - A stoichiometric amount of 0.1N HCl was added dropwise with swirling and gentle heating to dissolve the BZ (base). The solution was diluted to volume with double-distilled water, sealed in a serum vial, and stored under refrigeration. The BZ·HCl standards used for spiking the tap water were prepared by diluting the stock solution with double-distilled water.

#### 3.2.2 Preparation of Michler's Ketone Standards.

The stock solutions were prepared by dissolving Michler's Ketone in methyl ethyl ketone and sealing in serum vials using Teflon® lined septa. The vials were made opaque by wrapping in aluminum foil. The stock solutions were stored under refrigeration.

The Michler's Ketone standards that were used for the analysis by GC were prepared by dilution of the stock solution with methyl ethyl ketone. The standards were sealed in serum vials using Teflon® lined septa.

#### 3.2.3 Extraction of BZ from Tap Water.

The first extraction study for tap water determined the efficiency of the procedure for extracting BZ. One hundred grams of sodium sulfate were added to 1 liter of tap water in a 2-liter separatory funnel and dissolved by shaking for approximately 15 seconds. The sodium sulfate was used to increase the ionic strength of the water solution in an effort to "salt out" the BZ.

It was found that the amount of active chlorine (1.0 - 1.5 ppm) in the tap water, when made alkaline to pH 10-11, was sufficient to destroy BZ at the concentrations under study.

The active chlorine was removed by adjusting the pH of the sample to 3.0 using 10 ml of 1N  $\text{H}_2\text{SO}_4$  and followed by the addition of 1 ml of 0.2 percent sodium metabisulfite. The removal required less than 30 seconds for complete reaction.

After the chlorine level was checked using the test kit, BZ was added to the tap water in concentrations ranging from 1 ppb to 40 ppb.

To extract the BZ, the tap water was made basic (pH 10.5) by adding 100 ml of a 0.5 M sodium carbonate/0.1 M sodium hydroxide solution. Fifty milliliters of chloroform were added, and the mixture was shaken for 3 minutes. The layers were allowed to separate for 5 minutes before the chloroform layer was drawn off into a 250-ml separatory funnel. Another 50 ml of chloroform was added to the water in the 2-liter separatory funnel, and the mixture was again shaken for 3 minutes. The layers were allowed to separate for 10 minutes before the chloroform layer was drawn off and combined with the first chloroform extract.

Twenty-five milliliters of 0.001N  $\text{H}_2\text{SO}_4$  were added to the  $\text{CHCl}_3$  solution in the 250-ml separatory funnel, and the mixture was shaken for 3 minutes. The layers were allowed to separate for 5 minutes. The aqueous layer was transferred with a disposable Pasteur pipette to a 125-ml separatory funnel. Another 25 ml of 0.001N  $\text{H}_2\text{SO}_4$  were added to the  $\text{CHCl}_3$ , and the mixture was shaken for 3 minutes. The layers were allowed to separate for 5 minutes before the  $\text{CHCl}_3$  layer was drawn off and discarded. The aqueous layer was combined with the first aqueous extract.

The pH of the aqueous extract was adjusted to 11.0 using 1 ml of 0.5 M sodium carbonate. Five milliliters of  $\text{CHCl}_3$  were added, and the mixture was shaken for 3 minutes. The layers were allowed to separate for 5 minutes. The  $\text{CHCl}_3$  layer was drawn off into a conical centrifuge tube. Another 5 ml of  $\text{CHCl}_3$  layers were allowed to separate for 5 minutes before the  $\text{CHCl}_3$  layer was combined with the first extract.

The combined  $\text{CHCl}_3$  extracts were concentrated to approximately 0.1 ml in a constant-temperature water bath set at  $80^\circ\text{C}$  while purging with nitrogen. The nitrogen flow was controlled by a reflux timer connected to a solenoid. This alternatively started and stopped the nitrogen flow at 30-second intervals. The purpose of the interrupted flow was to set up a reflux to help wash down the walls of the centrifuge tube. When the concentrates reached a level of approximately 0.1 ml, they were removed from the bath, the internal standard (Michler's Ketone) was added, and the solutions were diluted to 0.2 ml with  $\text{CHCl}_3$ . The concentrates were stirred on a Vortex mixer before analyzing by GC.




Figure 1. Chromatogram of BZ Standard With Internal Standard




Figure 2. Chromatogram of Sample Extract  
Containing Contaminant Peak



Tables 1 and 2 present the BZ recovery from unaerated and aerated filters at the two BZ levels. The recovery shows the amount of BZ extracted from the filter as a percentage of the amount deposited on the filter. Both the peak height measurement and the area measurement recoveries are shown for the unaerated and aerated filters. The average recovery ( $\bar{X}$ ) and the standard deviation (S) are given for each measurement method. All recoveries were calculated by comparison to the peak height and areas of the external standard run before the sample.

Tables 3 and 4 present the recovery of the two levels of BZ from solution after concentration. This was done to check for loss of BZ during the concentration step. Recoveries were calculated as for tables 1 and 2.

#### 4.2 Water Studies.

Table 5 is a step-by-step outline of the method for water developed by CSL. An outline of the BCL method for air is also given for comparison. As can be seen, both methods are very similar and are complex, time consuming, and labor intensive. The prolonged procedure and long turnaround time makes the method not entirely suitable for most of the routine analyses in the BZ demilitarization plant.

Figures 3 and 4 illustrate the effect of chlorine on the recovery of BZ from tap water. Figure 3 is a chromatogram of an extraction of BZ from tap water without chlorine. Figure 4 shows the extraction of BZ from tap water containing chlorine. The BZ concentration in figure 4 is the same as in figure 3. In figure 3, BZ has a retention time of 3.10 minutes and Michler's Ketone has a retention time of 7.26 minutes. In figure 4, no peak appears at the retention time of BZ, but two extraneous peaks are evident at the retention times of 4.96 and 5.74 minutes. The peak of 7.30 minutes is Michler's Ketone. It was not determined whether the extraneous peaks were BZ degradation products.

Table 6 presents the recovery of BZ from tap water using the procedure described in 3.2.3. The recovery was calculated by comparison of the amount of BZ in the extract to a prepared GC standard. The GC standard had a BZ concentration that would approximate that in the extract if 100 percent of the BZ were recovered.

Depicted in table 7 is the recovery of BZ from tap water in the presence of contaminants, using the procedure described in 3.2.4. The effect on the recovery with and without the cyclohexane is shown. The recovery was calculated in the same manner as for table 6.

Table 8 shows the effect of salting on the recovery of BZ at the drinking water level (4 ppb). The recovery was calculated in the same manner as for table 6.

Tables 9 and 10 show the detection limit as calculated by the method of Hubaux and Vos.<sup>11</sup>

Figures 5 and 6 show the 90 percent confidence limits as calculated by the method of Hubaux and Vos.

Blank

Table 1. 82 Recovery at the 0.1 µg Level

Unaltered (N=25)				Acetated				1000 Liters (N=25)			
Sample	Peak Height (min)	Area (counts)	Recovery (%)	Sample	Peak Height (min)	Area (counts)	Recovery (%)	Sample	Peak Height (min)	Area (counts)	Recovery (%)
External Standard	92.5	1034	99.8	External Standard	82.0	1135	99.8	External Standard	78.0	1775	99.8
1	80.0	1032	86.5	1	49.0	566	49.9	1	34.5	1110	49.2
2	92.0	919	86.5	2	62.5	730	76.2	2	34.5	1006	44.2
External Standard	92.0	1003	99.8	External Standard	74.0	868	71.8	External Standard	73.5	1751	99.8
3	53.5	655	58.2	3	43.5	609	58.8	3	30.0	1065	60.8
4	66.0	778	71.7	4	49.0	652	66.2	4	31.0	1005	61.1
External Standard	96.0	806	99.8	External Standard	75.5	869	76.9	External Standard	71.5	1685	99.8
5	43.0	415	44.8	5	41.0	525	54.3	5	39.5	1412	55.2
6	28.0	297	28.2	6	32.0	306	42.4	6	46.0	1299	64.3
External Standard	106.0	905	99.8	External Standard	67.0	743	63.4	External Standard	114.5	2717	99.8
7	39.5	385	37.3	7	42.5	427	66.4	7	44.5	1041	38.9
8	56.0	595	54.7	8	44.5	454	61.1	8	49.0	1218	44.8
External Standard	98.0	915	99.8	External Standard	43.0	433	67.2	External Standard	119.0	2703	99.8
9	41.5	388	42.4	9	73.5	810	54.3	9	43.5	1777	53.4
10	49.0	432	50.0	10	43.0	440	58.5	10	44.0	1627	53.8
External Standard	97.0	1026	99.8	External Standard	56.0	545	76.2	External Standard	105.0	2827	99.8
11	63.0	693	63.0	11	51.5	517	70.1	11	51.0	1626	53.3
12	54.0	505	53.7	12	69.0	740	60.9	12	51.0	1497	53.0
External Standard	90.0	975	99.8	External Standard	42.0	482	62.3	13	48.5	1939	68.7
13	55.0	640	61.1	13	56.0	561	61.2	14	71.0	2021	71.6
14	61.5	690	68.3	14	43.0	410	55.4	15	72.5	1966	72.5
15	56.5	596	62.8	15	96.5	1124	47.1	16	71.5	1868	76.5
16	79.0	894	87.8	16	37.0	323	33.2	17	73.0	1995	73.2
External Standard	98.0	1106	99.8	External Standard	100.0	1095	42.5	18	73.0	1996	73.2
17	52.0	652	53.1	17	42.5	308	38.3	19	73.0	2322	73.2
18	66.5	779	67.9	18	121.5	1342	42.5	20	73.0	1634	61.2
External Standard	83.0	877	99.8	External Standard	39.0	340	32.1	21	73.0	1634	61.2
19	30.5	396	36.8	19	49.5	445	40.7	22	73.0	1634	61.2
20	33.0	363	42.2	20	32.5	473	43.2	23	73.0	1634	61.2
21	34.0	348	41.0	21	101.5	1169	43.2	24	73.0	1634	61.2
External Standard	86.5	821	99.8	External Standard	68.0	744	67.0	25	73.0	1634	61.2
22	26.5	358	25.9	22	43.6	426	37.4	26	73.0	1634	61.2
23	31.0	541	59.0	23	59.9	445	51.2	27	73.0	1634	61.2
External Standard	102.0	1058	99.8	External Standard	93.5	1027	93.3	28	73.0	1634	61.2
24	68.0	987	67.7	24	60.8	562	51.2	29	73.0	1634	61.2
25	62.0	733	60.8	25	69.3	562	60.8	30	73.0	1634	61.2
Σ (mean)				Σ (mean)				Σ (mean)			
s (standard deviation)				s (standard deviation)				s (standard deviation)			

Blank



Table 2. BZ Recovery at the 1.4 µg Level

Unreacted (N-25)					Aerated					1000 Liters (N-25)					Seawater				
Sample	Peak Height (mm)	Area (counts)	Peak Height	Area	Sample	Peak Height (mm)	Area (counts)	Peak Height	Area	Sample	Peak Height (mm)	Area (counts)	Peak Height	Area	Sample	Peak Height (mm)	Area (counts)	Peak Height	Area
External Standard	124.0	19540	-	-	External Standard	141.0	22490	-	-	External Standard	141.0	23280	-	-	External Standard	141.0	23280	-	-
1	28.5	13790	63.3	65.5	1	75.0	12170	53.2	54.1	1	75.0	12170	53.2	54.1	1	75.0	12170	53.2	54.1
2	81.0	14320	66.1	73.3	2	79.5	13980	56.4	62.2	2	79.5	13980	56.4	62.2	2	79.5	13980	56.4	62.2
External Standard	106.0	16410	-	-	External Standard	134.5	21670	-	-	External Standard	134.5	21670	-	-	External Standard	134.5	21670	-	-
3	87.0	13940	82.1	84.9	3	82.0	13510	61.0	62.3	3	82.0	13510	61.0	62.3	3	82.0	13510	61.0	62.3
4	74.0	12220	69.8	74.5	4	72.5	16140	53.9	74.5	4	72.5	16140	53.9	74.5	4	72.5	16140	53.9	74.5
External Standard	131.0	20860	-	-	External Standard	84.0	13610	62.4	62.8	External Standard	84.0	13610	62.4	62.8	External Standard	84.0	13610	62.4	62.8
5	75.0	12560	57.2	60.2	5	100.0	21340	63.9	64.2	5	100.0	21340	63.9	64.2	5	100.0	21340	63.9	64.2
6	76.0	11920	58.0	57.1	6	83.0	13700	69.6	70.3	6	83.0	13700	69.6	70.3	6	83.0	13700	69.6	70.3
7	82.5	16730	63.0	80.2	7	90.5	15000	76.2	76.9	7	90.5	15000	76.2	76.9	7	90.5	15000	76.2	76.9
External Standard	126.0	20270	-	-	External Standard	99.0	16420	75.0	78.1	External Standard	99.0	16420	75.0	78.1	External Standard	99.0	16420	75.0	78.1
8	91.5	18600	71.5	91.8	8	97.5	16660	75.0	78.1	8	97.5	16660	75.0	78.1	8	97.5	16660	75.0	78.1
9	106.5	22310	83.2	115	9	130.5	6208	-	-	9	130.5	6208	-	-	9	130.5	6208	-	-
External Standard	143.5	22860	-	-	External Standard	104.0	4914	79.7	79.2	External Standard	104.0	4914	79.7	79.2	External Standard	104.0	4914	79.7	79.2
10	81.0	13100	56.5	57.3	10	104.0	4914	79.7	79.2	10	104.0	4914	79.7	79.2	10	104.0	4914	79.7	79.2
External Standard	99.5	9028	-	-	External Standard	102.5	4821	78.5	77.7	External Standard	102.5	4821	78.5	77.7	External Standard	102.5	4821	78.5	77.7
11	56.0	3402	38.3	39.9	11	93.0	4366	71.3	70.3	11	93.0	4366	71.3	70.3	11	93.0	4366	71.3	70.3
12	77.0	7246	77.4	80.3	12	125.0	6140	93.0	94.6	12	125.0	6140	93.0	94.6	12	125.0	6140	93.0	94.6
13	62.5	5922	62.8	65.6	13	120.0	5807	75.2	76.0	13	120.0	5807	75.2	76.0	13	120.0	5807	75.2	76.0
External Standard	97.0	9084	-	-	External Standard	97.0	4666	60.5	55.6	External Standard	97.0	4666	60.5	55.6	External Standard	97.0	4666	60.5	55.6
14	57.0	4910	38.8	34.0	14	86.0	3412	68.2	68.1	14	86.0	3412	68.2	68.1	14	86.0	3412	68.2	68.1
15	61.0	5599	62.9	61.6	15	129.5	6178	78.8	78.4	15	129.5	6178	78.8	78.4	15	129.5	6178	78.8	78.4
16	65.5	6163	67.5	67.8	16	102.0	4846	83.8	82.8	16	102.0	4846	83.8	82.8	16	102.0	4846	83.8	82.8
17	51.0	4536	52.6	49.9	17	108.5	5114	70.3	69.7	17	108.5	5114	70.3	69.7	17	108.5	5114	70.3	69.7
External Standard	106.5	9680	-	-	External Standard	91.0	4305	-	-	External Standard	91.0	4305	-	-	External Standard	91.0	4305	-	-
18	48.5	4161	45.5	43.0	18	132.0	6263	79.2	80.9	18	132.0	6263	79.2	80.9	18	132.0	6263	79.2	80.9
19	53.5	4742	50.2	49.0	19	104.5	5145	79.2	79.4	19	104.5	5145	79.2	79.4	19	104.5	5145	79.2	79.4
20	62.0	5874	58.2	60.7	20	123.5	6051	74.1	73.5	20	123.5	6051	74.1	73.5	20	123.5	6051	74.1	73.5
21	58.5	5352	54.3	55.3	21	91.5	4445	91.5	90.7	21	91.5	4445	91.5	90.7	21	91.5	4445	91.5	90.7
External Standard	78.0	7166	-	-	External Standard	113.0	5489	81.8	79.8	External Standard	113.0	5489	81.8	79.8	External Standard	113.0	5489	81.8	79.8
22	64.0	6008	62.1	63.8	22	101.0	4289	-	-	22	101.0	4289	-	-	22	101.0	4289	-	-
23	34.0	5006	69.2	69.9	23	101.0	4289	-	-	23	101.0	4289	-	-	23	101.0	4289	-	-
24	97.0	4310	74.3	72.7	24	101.0	4289	-	-	24	101.0	4289	-	-	24	101.0	4289	-	-
25	92.0	4434	70.5	71.4	25	101.0	4289	-	-	25	101.0	4289	-	-	25	101.0	4289	-	-
Σ (mean)					Σ (mean)					Σ (mean)					Σ (mean)				
σ (standard deviation)					σ (standard deviation)					σ (standard deviation)					σ (standard deviation)				

Blank

Table 3. BZ Recovery After Concentration - 0.1  $\mu$ g Level

Sample	Peak Height (mm)	Area (counts)	Recovery (%)	
			Peak Height	Area
External Standard	98.0	5538	-	-
1	102.5	5379	105.0	97.1
2	89.5	4752	91.3	85.8
3	88.0	4592	89.8	82.9
4	79.5	4200	81.1	75.8
External Standard	94.5	5148	-	-
5	91.0	4837	96.3	94.0
6	101.0	5319	107.0	103.0
7	99.5	5205	105.0	101.0
8	83.0	4370	107.0	103.0
External Standard	105.5	5731	-	-
9	104.0	5459	98.6	95.3
10	97.0	5120	91.9	89.9
$\bar{x}$			95.3	91.0
s			8.57	8.63

Table 4. BZ Recovery After Concentration 1.4  $\mu$ g Level

Sample	Peak Height (mm)	Area (counts)	Recovery (%)	
			Peak Height	Area
External Standard	111.5	93500	-	-
	89.5	72160	80.3	77.2
	93.0	74460	83.4	79.6
	95.0	76180	85.2	81.5
	118.5	95160	106.0	102.0
External Standard	104.5	86740	-	-
	90.0	72560	86.1	83.6
	109.0	87520	104.0	101.0
External Standard	89.5	74180	-	-
	97.0	79080	108.0	107.0
	92.0	74640	103.0	101.0
	111.0	90780	124.0	122.0
	106.0	86380	118.0	116.0
$\bar{x}$			99.8	97.1
s			15.3	15.8

Table 5. Comparison of Methods

Step #	Method Developed by CSL		Method Developed by BCL	
	Operation	Estimated Time (min)	Operation	Estimated Time (min)
1	Add H <sub>2</sub> SO <sub>4</sub> and dissolve	1	Spike filter with BZ-d <sub>10</sub>	1
2	Adjust pH to 3.0	0.25	Add 0.1N HCl or acetic acid	0.5
3	Add 0.2% Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	0.25	Shake	2
4	Let solution stand	0.5	Add cyclohexane	0.5
5	Check chlorine level	1	Shake	2
6	Add cyclohexane	0.5	Let settle	2
7	Shake	1	Remove cyclohexane	0.5
8	Allow layers to separate	5	Add cyclohexane	0.5
9	Draw off water layer into 2-liter separatory funnel	0.5	Shake	2
10	Adjust pH to 10.5	0.25	Let settle	2
11	Add CHCl <sub>3</sub>	0.5	Remove cyclohexane	0.5
12	Shake	3	Add CHCl <sub>3</sub>	0.5
13	Allow layers to separate	5	Add Na <sub>2</sub> CO <sub>3</sub> to pH 11 (check pH)	2
14	Draw off CHCl <sub>3</sub> layer into 250ml separatory funnel	0.5	Shake	2
15	Add fresh CHCl <sub>3</sub>	0.5	Filter into separatory funnel	2
16	Shake	3	Let settle	2
17	Allow layers to separate	10	Remove CHCl <sub>3</sub> into 250ml Erlenmeyer	0.5
18	Draw off CHCl <sub>3</sub> into 250ml separatory funnel	0.5	Add fresh CHCl <sub>3</sub>	0.5
19	Add 0.001N H <sub>2</sub> SO <sub>4</sub>	0.5	Shake	2
20	Shake	3	Let settle	2
21	Allow layers to separate	5	Remove CHCl <sub>3</sub> into same Erlenmeyer	0.5
22	Transfer H <sub>2</sub> SO <sub>4</sub> to 250ml separatory funnel	2	Add fresh CHCl <sub>3</sub>	0.5
23	Add fresh 0.001N H <sub>2</sub> SO <sub>4</sub>	0.5	Shake	2
24	Shake	3	Let settle	2
25	Allow layers to separate	5	Remove CHCl <sub>3</sub> into same Erlenmeyer	0.5
26	Draw off & discard CHCl <sub>3</sub> layer	0.5	Add K <sub>2</sub> CO <sub>3</sub> (to dry)	0.5
27	Draw off H <sub>2</sub> SO <sub>4</sub> into 125ml separatory funnel	0.5	Filter into 250ml rb flask	2
28	Adjust pH to 11.0	0.25	Rinse (2x)	1
29	Add CHCl <sub>3</sub>	0.25	Attach to rotovaporator	1
30	Shake	3	Evaporate	15
31	Allow layers to separate	5	Remove from rotovaporator	1
32	Draw off CHCl <sub>3</sub> into centrifuge tube	0.5	Pour into K <sub>2</sub> O concentrator	1
33	Add fresh CHCl <sub>3</sub>	0.25	Rinse (2x)	1
34	Shake	3	Concentrate	15
35	Allow layers to separate	5	Transfer to vial	0.5
36	Draw off CHCl <sub>3</sub> into centrifuge tube	0.5	Seal vial	0.5
37	Evaporate CHCl <sub>3</sub> to 0.1ml	15	Place in autosampler rack	0.5
38	Add internal standard	0.25		
39	Dilute to 0.2ml	0.25		
	TOTAL	86.5	TOTAL	71.5




Figure 3. Chromatogram of Extraction of BZ from Tap Water




Figure 4. Chromatogram Showing Effect of Chlorine on Extraction of BZ from Tap Water

Blank

Table 4. Recovery of BZ from Tap Water

[illegible]

Blank



Table 7. Recovery of BZ From Tap Water in the Presence of Contaminants

Without Cyclohexane Wash						
BZ Concentration (ppb)	Contaminant	Contaminant Concentration (ppm)	Recovery (%)	$\bar{X}$	S	%S
1.24	Sodium Benzilate	11.2	0 0	---	---	---
1.24	Benzophenone	10.0	0 0	---	---	---
1.24	3-Quinuclidinol	7.75	45.4 40.0 41.7	42.4	2.76	6.51
With Cyclohexane Wash						
BZ Concentration (ppb)	Contaminant	Contaminant Concentration (ppm)	Recovery (%)	$\bar{X}$	S	%S
1.24	---	---	44.7 37.3 29.3	37.1	7.70	20.8
1.24	Sodium Benzilate	11.2	16.6 29.2 21.1	22.3	6.39	28.7
1.24	Sodium Benzilate	1.12	33.6 21.4 19.1	24.7	7.79	31.5
1.24	Benzophenone	10.0	13.8 34.6 16.3	21.6	11.4	52.8
1.24	Benzophenone	1.00	17.5 24.8 25.4	23.8	4.39	18.4

Table 8. Effect of Salting on Recovery of BZ

<u>Without Salting</u>			
<u>BZ Concentration (ppb)</u>	<u>Michler's Ketone Peak Height (mm)</u>	<u>BZ Peak Height (mm)</u>	<u>Recovery (%)</u>
Standard (3.98)	67.5	71.0	-
4.0	63.5	22.5	33.7
4.0	79.5	24.0	28.5
Standard (3.98)	66.0	72.5	-
4.0	81.0	33.0	37.0
		-	
		x	33.1
		s	4.29
<u>With Salting</u>			
<u>BZ Concentration (ppb)</u>	<u>Michler's Ketone Peak Height (mm)</u>	<u>BZ Peak Height (mm)</u>	<u>Recovery (%)</u>
Standard (4.02)	69.5	72.0	-
4.0	60.0	32.75	53.0
4.0	80.5	46.25	55.8
Standard (4.02)	67.5	74.75	-
4.0	54.5	29.0	48.4
Standard (3.98)	66.0	72.5	-
4.0	68.0	38.5	51.5
		-	
		x	52.2
		s	3.08

Table 9. Hubaux - Vos Detection Limit

Target Value (ppb)	Found Value (ppb)
40	34.6
40	32.5
40	32.3
40	32.7
40	32.2
40	31.8
40	31.0
40	33.6
40	33.7
40	34.4
4	2.12
4	2.23
4	1.94
4	2.06
4	2.10
4	1.96
4	2.08
2	0.77
2	0.78
2	0.70
2	0.73
2	0.71
1	0.48
1	0.37
1	0.36
1	0.52
1	0.46
1	0.34
0	0
0	0
0	0
0	0

Detection Limit  $x(d) = 3.2$  ppb

Table 10. Hubaux-Vos Detection Limit

Target Value (ppb)	Found Value (ppb)
4	2.12
4	2.23
4	1.94
4	2.06
4	2.10
4	1.96
4	2.08
2	0.77
2	0.78
2	0.70
2	0.73
2	0.71
1	0.48
1	0.37
1	0.36
1	0.52
1	0.46
1	0.34
0	0
0	0
0	0
0	0

Detection Limit  $x(d) = 0.95$




Figure 5. 90 Percent Confidence Limit




Figure 6. 90 Percent Confidence Limit

## 5. DISCUSSION OF RESULTS

### 5.1 Aerosol Studies.

The BZ recoveries were lower and more variable than would be desired for routine laboratory analysis. The low and variable recovery from the filters could be due to the polarity of the glass fiber filters which bound the BZ to the filters and the inability of dioxane to remove it. Deiner, Herd, and Vigus, Munitions Division, CSL,\* showed that 100 percent recoveries from glass-fiber aerosol filters were attainable using an aqueous acidic medium. It is concluded that dioxane employed as a solvent is not sufficiently polar to afford complete removal of BZ from the filter.

In addition, much of the variability could be due to the NPD. The detector is very sensitive and selective for nitrogen and phosphorus compounds, but it can also be nonreproducible in its response to a sample. The NPD does not show long-term response stability. Attempts to quantitate BZ in a sample by comparing it to a previously injected standard may produce variable results. It is now thought that the lack of a suitable internal standard to help overcome the effects of the detector instability may be the cause for much of the variation in the results of the tests conducted. As a result of BCL's efforts, it was shown that Michler's Ketone could be successfully employed as an internal standard<sup>12</sup> with the NPD.

The results of the concentration tests did not show any appreciable loss of BZ at either dosage level.

#### 5.1.1 Statistical Analysis of the Data.

A statistical analysis of the data was done to determine if there was any significant difference at the two BZ levels between the amount of BZ recovered from the aerated filters and the amount of BZ recovered from the unaerated filters.

The data obtained from peak height measurements and from area measurements were analyzed separately.

##### 5.1.1.1 Peak Height Measurements.

The frequency distribution of the peak height measurement data from tables 1 and 2 are shown in figure 7.

#### Analysis of Variance.

Based upon Bartlett's test shown below, there would be a 17.2 percent risk of error if it were concluded that a significant difference did exist in the variability of the different samples. It was concluded that the variances of the samples were equivalent regardless of the aeration or BZ level used.

\*Deiner, A., Herd, R. E., and Vigus, E. S. The Utilization of an Ion-Pairing Reaction to Trace Concentrations of 3-Quinuclidinyl Benzilate (BZ) and Application to Demilitarization Operations. (unpublished data, June 1983).




Figure 7. FREQUENCY DIAGRAM FOR PEAK HEIGHT DATA OF TABLES 1 & 2



### BARTLETT'S TEST

$$\begin{aligned} M &= 2.3025 (\nu \log (\sum_i S_i^2 / \nu - \sum \nu_i \log S_i^2)) \\ &= 7.853 \end{aligned}$$

$$\begin{aligned} C &= 1 + \frac{1}{3(g-1)} \left( \frac{1}{\nu_i} - \frac{1}{\nu} \right) \\ &= 1.016 \end{aligned}$$

$$M/C = 7.728 = \chi^2_{0.172, 5}$$

$$\nu = \sum \nu_i$$

$$\nu_i = \text{degrees of freedom of the } i_{th} \text{ sample.}$$

$$S_i^2 = \text{independent estimate of the } i_{th} \text{ variance.}$$

$$g = \text{number of variances compared.}$$

$$\chi_{\alpha, g-1} = \text{percentage point of the Chi-Square distribution at the } \alpha \text{ point and } g-1 \text{ degrees of freedom.}$$

Assuming a completely cross-classified, random effects model, the analysis of the variance presented in table 11 shows that:

- (1) There would be a 31.5 percent risk of error if it were that aeration had an effect on the process.
- (2) There would be less than 0.05 percent risk of error if it were concluded that the BZ level had an effect on the process.
- (3) There would be a 9.5 percent risk of error if it were concluded that there was an effect on the process due to the interaction of the BZ level and the aeration volume.

From (1) above, it was concluded that aeration may have an effect on the amount of BZ recovered.

From (2) above, it was concluded that the level of BZ had no effect on the percent of recovery.

From (3) above, it was decided that more testing would be necessary before any firm conclusions could be drawn concerning the effect of the BZ level or the effect of the interaction of the BZ level and the aeration volume on the amount of BZ recovered.

Table 11. Analysis of Variance of Peak Height Measurement Data

Source	Degree of freedom	Sum of squares	Mean square	Mean square ratio	Level for $H_0$ rejection
BZ Level	1	3247.096	3247.096	20.922	<0.0005
Aeration	2	358.0769	179.038	1.154	0.315
Interaction	2	768.9266	364.463	2.348	0.095
Residual	144	22347.8088	155.193		
Total	149	26681.9083			

5.1.1.2 Sample Size.

The residual mean square shown in table 11 is an estimate of the experimental error of the peak height measurement method. Assuming that the residual mean square is the population variance ( $S^2$ ), the number of GC injections (n) needed to reliably estimate the population can be calculated using the following equation:

$$n = (t_{1-\alpha/2, \phi}^2 \hat{S}^2) / d^2$$

$t_{1-\alpha/2, \phi}$  -  $(1-\alpha/2)$  percentile of the t distribution at  $\phi$  degrees of freedom ( $\phi = 144$ )

$\hat{S}^2$  - estimate of the population variance ( $\hat{S}^2 = 155.193$ )

d - error in the estimate of the population mean

Table 12 shows the calculated number of GC injections needed depending upon the 100 (1- $\alpha$ ) percent confidence level and the interval  $\bar{X} \pm d$  percent one is willing to accept. It can be seen from the table that if the average value for a sample is to be within  $\pm 10$  percent of the true average value ( $\bar{X}$ ) 95 percent of the time, then at least six injections would be necessary.

Table 12. Sample Size

100 (1- $\alpha$ )% Confidence	n		
	d = 5	d = 10	d = 15
99	43	11	5
98	35	9	4
95	25	6	3
90	17	5	2
80	11	3	2

### 5.1.1.3 Area Measurement.

The frequency distributions of the area measurement data from tables 1 and 2 are shown in figure 8.

An analysis of variance was not done on the area measurement data because the assumption of homogeneity of variance between samples was proven not valid by a Bartlett's Test for the homogeneity of variances (shown below). The results of this test showed that there would be a 1.1 percent risk of error if it were concluded that a significant difference exists in the variability of the different samples. It was, therefore, concluded that the variability between samples was not constant.

#### BARTLETT'S TEST

$$M = 2.3026 (\sum v_i \log (\sum v_i S_i^2 v) - \sum v_i \log s_i^2) \\ = 15.142$$


$$C = 1 + \frac{1}{3(g-1)} \left( \frac{1}{v_i} - \frac{1}{v} \right) \\ = 1.016$$

$$M/C = 14.901 = \chi^2_{0.11,5}$$

### 5.1.1.4 t-Test and F-Test.

Differences in the mean and variance between individual samples were checked by t-tests and F-tests (see tables 13 and 14). While differences were found to exist between the mean and variances of individual samples, no pattern was found which would indicate an effect on the amount of BZ recovered by either the aeration volume or the BZ level used.

Table 13. t-Test of Area Measurement Data







Figure 8. FREQUENCY DIAGRAM FOR AREA DATA OF TABLES 1 & 2

Table 14. F-Test of Area Measurement Data



Area Measurement Method Versus Peak Height Measurement Method. Statistical analysis of the data showed that the peak height measurement method was more reliable than the area measurement method. The authors believe that this is true because it was possible to visually determine where the BZ peak began, but the change in slope frequently was not sufficient to trigger the electronic integrator at the true beginning and end of the BZ peak. This caused erroneous area count readings.

## 5.2 Water Studies.

The presence of active chlorine in water presents a potential problem for accurate analysis. Though BZ is stable to oxidants at pH 3-4<sup>13</sup> used in this assay, at a pH greater than 7, BZ at the 10-fold drinking water level is rapidly degraded. This is aptly illustrated in figures 3 and 4. The procedure developed to measure and destroy active chlorine is rapid and sufficiently accurate to preclude degradation of BZ samples prior to analysis.

Shown in figure 9 is the average recovery of BZ from tap water at various spike levels. The low recoveries at the low spike levels are due to intrinsic solubility of BZ in water.

From the data base,<sup>7</sup> the solubility of BZ base in water ranges from 11.8 ppm to 540 ppm. The data base also cites the solubility of BZ in  $CCl_4$  as 4.5 gm/liter while the solubility in  $CHCl_3$  is cited as 130 gm/liter. The distribution coefficient listed for  $CCl_4$  is approximately  $3.3 \times 10^2$ . One may expect a far more favorable distribution coefficient with  $CHCl_3$  due to the greater solubility of BZ in this solvent. From the data at the lowest level, the total loss of BZ is approximately 0.58 ppb. This reflects the ultimate efficiency one may expect using the procedure.

Blank

Blank

The losses observed at the low levels of spiking were minimized by "salting" the BZ out of aqueous solution using 10 percent  $\text{Na}_2\text{SO}_4$ . Shown in table 8 are the results of using  $\text{Na}_2\text{SO}_4$  to "salt out" BZ. This effect at the low level is appreciable and results in being able to quantitate BZ at one-quarter of the drinking water level (1ppb).

It was anticipated that under field conditions, large levels of contaminants would be present in process water. As part of this work, the effects of high levels of contaminants present along with the BZ were determined using the cleanup method in section 3.2.4.

The three contaminants studied were 3-quinuclidinol, sodium benzilate, and benzophenone. Sodium benzilate is highly soluble in water. Under acidic conditions such as in scrubbing for  $\text{SO}_2$  and in venturi washing, it is converted to benzilic acid. Benzilic acid has a solubility of 168.2 mg/100 ml (1682 ppm).<sup>\*</sup> This high solubility could result in benzilic acid being present in the final extract.

The anticipated degradation product of BZ due to incineration is benzophenone. This material has a solubility of 6.9 mg/100 ml (690 ppm).<sup>\*</sup> Though benzophenone is a neutral, nonionic species, its inherent high solubility in water could result in some benzophenone not being completely removed by the initial cyclohexane wash.

Early in our study we found that benzophenone and benzilic acid, injected on-column with BZ, results in almost complete degradation of the BZ. Monsanto, under contract DA18-035-AMC-136(A)<sup>14</sup> and Sass<sup>15</sup> showed that the rate of degradation of BZ was increased when benzophenone or benzilate was present. BCL has shown that benzilate degrades to benzophenone when injected on-column.<sup>14</sup>

A study was performed to measure the adequacy of the cleanup technique to remove this interference. Table 7 shows that removal of these interferences by the extraction technique is reasonably adequate. Some losses of BZ were observed, but total destruction was prevented.

The detection limit for tap water, as determined by the method of Hubaux and Vos,<sup>15</sup> was found to be 3.2 ppb when all data are included (table 9). When the values at 40 ppb are excluded, a detection limit of 0.95 ppb is obtained (table 10). Detection limits were calculated at the 90 percent confidence level.

## 6. CONCLUSIONS

### 6.1 Employment of NPD for BZ Analysis.

The manufacturer of the NPD recommends that chlorinated solvents not be used. BZ can be quantitatively removed from aerosol filters with acidic aqueous solvents. Extraction of the BZ from the aqueous phase must be performed with neutral, low-boiling, halogen-free, water-immiscible, organic solvents if the NPD is employed. The use of alcohols, aldehydes, ketones, and esters is precluded due to their water solubility; leaving neutral hydrocarbons as the only viable extraction



solvents. Unfortunately, BZ has a low solubility in neutral hydrocarbons. Due to the intrinsic solubility of BZ in water, it is believed that even with multiple extractions of the aqueous phase with hydrocarbon, extraction efficiencies approaching those obtained with a chlorinated hydrocarbon would be hard to achieve. In addition, the inherent instability of the NPD precludes its use on a routine basis.

It is recommended that the NPD not be employed for routine analysis of BZ.

## 6.2 Employment of FID for BZ Analysis.

Our studies show that sodium metabisulfite is an effective reducing agent to remove active chlorine from process water. In order to preserve water samples, it is recommended that the samples be acidified to a pH of 4, and that the active chlorine be removed by the process described. At this pH, adsorption of BZ on the glass walls of containers is negligible. As a further precaution, the samples should be stored at 10°C or lower if there are long delays in analysis.

It has been the experience of these laboratory personnel that the routine use of an internal standard to quantitate trace levels of BZ by GC with any detector is a fundamental requirement. Day-to-day variations in the GC detector sensitivity require the employment of a standard directly in the solution being analyzed. We find that Michler's Ketone is a suitable material as a reference standard. It is easy to purify, highly stable in solution, and does not react with BZ. The only drawback is that the retention time of Michler's Ketone is longer than desired. However, this drawback does not appreciably add to the length of the analysis.

The effect of interferences on the analysis was investigated. The cleanup technique effectively removes these materials and allows the quantification of BZ to one-quarter the drinking water level (1 ppb).

The procedures developed at CSL and at BCL are virtually identical. The method must be performed by highly-trained personnel who understand the necessity of carrying out the complex steps in the analytical scheme correctly and diligently. The procedure is long and labor-intensive. It appears that there are no means of shortening the analysis. Carried out correctly, the method is capable of quantitating BZ in water with good precision. It is recommended that the procedures developed at CSL and BCL (using GC-MS with BZd<sub>10</sub> as the internal standard) be employed to analyze process water, brines, and scrubber effluents. Presently, we do not feel that GC analysis is the best analytical tool to analyze a large number of filter samples. It is recommended that another analytical method be sought to analyze the large number of filter samples that will be generated every day at the actual demilitarization site.

## LITERATURE CITED

1. Skinner, V. L. Concept Plan for Demilitarization of Incapacitating BZ Agent/Munition, Office of the Project Manager for Chemical Demilitarization and Installation Restoration (presently USATHAMA), Aberdeen Proving Ground, Maryland. October 1976.
2. Sass, S., Pinsky, S., Schlotzhauer, W., and Beitsch, N. Chemical Research and Development Laboratories Report CRDLR 3183. Basic Esters of Glycolic Acids. VII. Gas-Liquid Chromatography of BZ and Process Intermediated. September 1963. UNCLASSIFIED Report.
3. Ellin, R. I. Chemical Research and Development Laboratories Technical Memorandum CRDL 23-27. The Use of Gas-Liquid Chromatography for the Quantitative Determination of BZ in Biological Systems. July 1962. UNCLASSIFIED Report.
4. Demek, M. M., and Epstein, J. Chemical Research and Development Laboratories Report CRDLR 3186. Chemistry of BZ. II. Estimation of BZ by Reaction with Tropaeolin 00. October 1963. UNCLASSIFIED Report.
5. Rosenblatt, D. H., Demek, M. M., and Epstein, J. Chemical Research and Development Laboratories Report CRDLR 3151. Chemistry of BZ. I. Reaction of BZ with Iodine in Aqueous and Organic Solution. November 1962. UNCLASSIFIED Report.
6. Petersen, B. A., Riggins, R. M., Shafer, K. H., Wyant, R. E., and Graffeo, A. P. Final Report, Contract DAAK40-73-C-0142. Evaluation of Analytical Methods for the Determination of BZ. December 1977. UNCLASSIFIED Report.
7. Rosenblatt, D. H., Dacre, J. C., Shiotsuka, R. N., and Rowlett, C. D. US Army Medical Bioengineering Research and Development Laboratory Technical Report 7710. Problem Definition Studies on Potential Environmental Pollutants. VII. Chemistry and Toxicology of BZ (3-Quinuclidinyl Benzilate). August 1977. UNCLASSIFIED Report.
8. Miller, B. A., and Junglaeus, G. A. Test Report - Task 2, Contract DAAK11-78-C-0096. Comparative Evaluation of Gas Chromatographic Detectors. November 1979. UNCLASSIFIED Report.
9. The Merck Index, 9th Edition, Merck & Co., Inc., Rahway, NJ. 1976.
10. Battelle Columbus Laboratories. Monthly Report, Contract DAAK11-78-C-0096. Monthly Progress Report on Engineering and Technical Support of Agent BZ Disposal Process. 22 August - 30 September 1978. UNCLASSIFIED Report.

11. Hubaux, A., and Vos, G. Decision and Detection Limits for Linear Calibration Curves. Anal. Chem. 42 (8), 849 (1970).

12. Battelle Columbus Laboratories. Monthly Report, Contract DAAK11-78-C-0096. Monthly Progress Report on Engineering and Technical Support of Agent BZ Disposal Process. 1-31 October 1978. UNCLASSIFIED Report.

13. Deiner, A., Kipp, R., and Herd, R. E. Chemical Research and Development Laboratories Technical Memorandum CRDL 33-38. Field Sampling and Analysis of EA2277. August 1961. UNCLASSIFIED Report.

14. Richardson, G. A. Final Report, Contract DA18-035-AMC-136(A). Physiocochemical Property-Structure Relationship of Glycolic Esters. November 1966. UNCLASSIFIED Report.

15. Sass, S., and Davis, R. M. Chemical Research and Development Laboratories Report CRDLR 3232. Basic Esters of Glycolic Acids. IX. Thermogravimetric and Other Thermal Studies on Some Glycolates. December 1964. UNCLASSIFIED Report.

# DISTRIBUTION LIST 3

Names	Copies	Names	Copies
CHEMICAL RESEARCH AND DEVELOPMENT CENTER		DEPARTMENT OF THE ARMY	
ATTN: DRSMC-CLB (A)	1		
ATTN: DRSMC-CLB-C (A)	1	HQDA	
ATTN: DRSMC-CLB-PO (A)	1	ATTN: DAMO-NCC	1
ATTN: DRSMC-CLB-R (A)	1	WASH DC 20310	
ATTN: DRSMC-CLB-R(A) (A)	1		
ATTN: DRSMC-CLB-R(S) (A)	1	Federal Emergency Management Agency	
ATTN: DRSMC-CLB-T (A)	1	Office of Research/NPP	
ATTN: DRSMC-CLB-TE (A)	1	ATTN: David W. Bensen	1
ATTN: DRSMC-CLF (A)	1	Washington, DC 20472	
ATTN: DRSMC-CLJ-R (A)	1		
ATTN: DRSMC-CLJ-L (A)	2	HQDA	
ATTN: DRSMC-CLJ-M (A)	1	Office of the Deputy Chief of Staff for	
ATTN: DRSMC-CLJ-P (A)	1	Research, Development & Acquisition	
ATTN: DRSMC-CLN (A)	1	ATTN: DAMA-CSS-C	1
ATTN: DRSMC-CLN-D (A)	1	Washington, DC 20310	
ATTN: DRSMC-CLN-TE (A)	1		
ATTN: DRSMC-CLT (A)	1	HQ Sixth US Army	
ATTN: DRSMC-CLW-C (A)	1	ATTN: AFKC-OP-NBC	1
ATTN: DRSMC-CLY-A (A)	1	Presidio of San Francisco, CA 94129	
ATTN: DRSMC-CLR-I (A)	1		
ATTN: DRSMC-CLY-R (A)	1	Commander	
		USASTCFEO	
COPIES FOR AUTHOR(S)		ATTN: MAJ Mikeworth	1
DRSMC-CLN-D (A)	12	APO San Francisco 96328	
RECORD COPY: DRSMC-CLN-A (A)	1		
DEPARTMENT OF DEFENSE		Commander	
		USA Nuclear & Chemical Agency	
Defense Technical Information Center		ATTN: MONA-WE	1
ATTN: DTIC-DDA-2	2	7500 Backlick Rd, Bldg 2073	
Cameron Station, Building 5		Springfield, VA 22150	
Alexandria, VA 22314			
Director		Army Research Office	
Defense Intelligence Agency		ATTN: DRXRO-CB (Dr. R. Ghirardelli)	1
ATTN: DB-4G1	1	P.O. Box 12211	
Washington, DC 20301		Research Triangle Park, NC 27709	
Commander			
USASED, USAINSCOM			
ATTN: IAFM-SED-III	1		
Fort Meade, MD 20755			

OFFICE OF THE SURGEON GENERAL

Commander

USA Medical Research Institute of  
Chemical Defense

ATTN: SGRD-UV-L

Aberdeen Proving Ground, MD 21010

Commander

US Army Environmental Hygiene Agency

ATTN: HSHB-O (B. Donovan)

Aberdeen Proving Ground, MD 21010

US ARMY MATERIEL DEVELOPMENT AND  
READINESS COMMAND

Commander

HQ, DARCOM

ATTN: DRCED (BG Robinson)

5001 Eisenhower Ave

Alexandria, VA 22333

Commander

USA Materiel Development and  
Readiness Command

ATTN: DRCSE-P

5001 Eisenhower Ave

Alexandria, VA 22333

Project Manager Smoke/Obscurants

ATTN: DRCPM-SMK-S

Aberdeen Proving Ground, MD 21005

Commander

USA Toxic & Hazardous Materials Agency

ATTN: DRXTH-ES

ATTN: DRXTH-TE

Aberdeen Proving Ground, MD 21010

Commander

USA Foreign Science & Technology Center

ATTN: DRXST-MT3

220 Seventh St., NE

Charlottesville, VA 22901

Director

USA Materiel Systems Analysis Activity

ATTN: DRXSY-MP

ATTN: DRXSY-CR (Mr. Metz)

Aberdeen Proving Ground, MD 21005

Commander

USA Missile Command

Redstone Scientific Information Center

ATTN: DRSMI-RPR (Documents)

Redstone Arsenal, AL 35898

Director

DARCOM Field Safety Activity

ATTN: DRXOS-C

Charlestown, IN 47111

Commander

USA Natick Research and Development  
Laboratories

ATTN: DRDNA-O

ATTN: DRDNA-IC

ATTN: DRDNA-IM

ATTN: DRDNA-ITF (Dr. Roy W. Roth)

Natick, MA 01760

US ARMY ARMAMENT, MUNITIONS AND  
CHEMICAL COMMAND

Commander

USA Armament, Munitions and Chemical Command

ATTN: DRSMC-ASN (R)

ATTN: DRSMC-IRW (R)

ATTN: DRSMC-LEP-L (R)

ATTN: DRSMC-SF (R)

Rock Island, IL 61299

Commander

USA AMCCOM

ATTN: SMCTE

Aberdeen Proving Ground, MD 21010

Commander

USA Dugway Proving Ground

ATTN: Technical Library (Docu Sect)

Dugway, UT 84022

US ARMY ARMAMENT RESEARCH AND  
DEVELOPMENT CENTER

Commander

USA Armament Research and  
Development Center

ATTN: DRSMC-LCA-L (D) 1  
ATTN: DRSMC-LCE-C (D) 1  
ATTN: DRSMC-LCU-CE (D) 1  
ATTN: DRSMC-SCA-C (D) 1  
ATTN: DRSMC-SCM (D) 1  
ATTN: DRSMC-SCP (D) 1  
ATTN: DRSMC-SCS (D) 1  
ATTN: DRSMC-TDC (D) (Dr. D. Gyorg) 1  
ATTN: DRSMC-TSS (D) 2  
ATTN: DRCPM-CAWS-AM (D) 1  
Dover, NJ 07801

USA Armament Research and  
Development Center

Resident Operations Office

ATTN: DRSMC-TSE-OA (Robert Thresher) 1  
National Space Technology Laboratories  
NSTL Station, MS 39529

Commander

ARRADCOM

ATTN: DRSMC-QAC-E 1  
Aberdeen Proving Ground, MD 21010

Director

USA Ballistic Research Laboratory

ARRADCOM

ATTN: DRSMC-TSB-S 1  
Aberdeen Proving Ground, MD 21005

Commander

USA Technical Detachment

US Naval EOD Technology Center

Indian Head, MD 20640

US ARMY TRAINING & DOCTRINE COMMAND

Commandant

USA Infantry School

ATTN: CTDD, CSD, NBC Branch 1  
Fort Benning, GA 31905

Commandant

USA Missile & Munitions Center  
and School

ATTN: ATSK-CM 1

ATTN: ATSK-TME 1

Redstone Arsenal, AL 35898

Commander

USA Logistics Center

ATTN: ATCL-MG 1

ATTN: DLSIE 1

Fort Lee, VA 23801

Commandant

USA Chemical School

ATTN: ATZN-CM-C 1

ATTN: ATZN-CM-AFL 2

ATTN: ATZN-CM-TPC 2

Fort McClellan, AL 36205

Commander

USA Infantry Center

ATTN: ATSH-CD-MS-C 1

Fort Benning, GA 31905

Commander

USA Infantry Center

Directorate of Plans & Training

ATTN: ATZB-DPT-PO-NBC 1

Fort Benning, GA 31905

Commander

USA Training and Doctrine Command

ATTN: ATCD-N 1

Fort Monroe, VA 23651

Commander

USA Armor Center

ATTN: ATZX-PPT-PO-C 1

Fort Knox, KY 40121

Commander

USA Combined Arms Center and  
Fort Leavenworth

ATTN: ATZL-CAM-IM 1

Fort Leavenworth, KS 66027

US ARMY TEST & EVALUATION COMMAND

Commander

USA Test & Evaluation Command

ATTN: DRSTE-CM-F

ATTN: DRSTE-CT-T

Aberdeen Proving Ground, MD 21005

DEPARTMENT OF THE NAVY

Project Manager

Theatre Nuclear Warfare Project Office

ATTN: PM-23 (Dr. Patton)

ATTN: TN-09C

Navy Department

Washington, DC 20360

Chief of Naval Research

ATTN: Code 441

800 N. Quincy Street

Arlington, VA 22217

Commander

Naval Explosive Ordnance Disposal

Technology Center

ATTN: AC-3

Indian Head, MD 20640

Officer-in-Charge

Marine Corps Detachment

Naval Explosive Ordnance Disposal

Technology Center

Indian Head, MD 20640

Commander

Naval Surface Weapons Center

Code G51

Dahlgren, VA 22448

Chief, Bureau of Medicine & Surgery

Department of the Navy

ATTN: MED 3C33

Washington, DC 20372

Commander

Naval Air Development Center

ATTN: Code 2012 (Dr. Robert Helmbold)

Warminster, PA 18974

US MARINE CORPS

Commanding General

Marine Corps Development and

Education Command

ATTN: Fire Power Division, D091

Quantico, VA 22134

DEPARTMENT OF THE AIR FORCE

Department of the Air Force

Headquarters Foreign Technology Division

ATTN: TQTR

Wright-Patterson AFB, OH 45433

ASD/AESD

Wright-Patterson AFB, OH 45433

AFAMRL/HE

ATTN: Dr. Clyde Reploggle

Wright-Patterson AFB, OH 45433

AFAMRL/TS

ATTN: COL Johnson

Wright-Patterson AFB, OH 45433

AFWAL/FIEEC (Wendell Banks)

Wright-Patterson AFB, OH 45433

HQ AFSC/SDZ

ATTN: CPT D. Riediger

Andrews AFB, MD 20334

HQ, AFSC/SGB

Andrews AFB, DC 20334

HQ, NORAD

ATTN: J-3TU

Peterson AFB, CO 80914

USAF TAWC/THL

Eglin AFB, FL 32542

AFATL/DLV

Eglin AFB, FL 32542



USAF SC  
ATTN: AD/YQ  
ATTN: AD/YQO (MAJ Owens)  
Eglin AFB, FL 32542

AD/XRO  
Eglin AFB, FL 32542

OUTSIDE AGENCIES

Battelle, Columbus Laboratories  
ATTN: TACTEC  
505 King Avenue  
Columbus, OH 43201

Toxicology Information Center, JH 652  
National Research Council  
2101 Constitution Ave., NW  
Washington, DC 20418

US Public Health Service  
Center for Disease Control  
ATTN: Lewis Webb, Jr.  
Building 4, Room 232  
Atlanta, GA 30333

Director  
1 Central Intelligence Agency  
1 ATTN: AMR/ORD/DD/S&T  
Washington, DC 20505

1 ADDITIONAL ADDRESSEES

Commandant  
Academy of Health Sciences, US Army  
ATTN: HSHA-CDH  
ATTN: HSHA-IPM  
1 Fort Sam Houston, TX 78234  
2